



# Transforming the Landscape of Rehabilitation





### Transforming the Landscape of Rehabilitation

At MedsEasy, we deeply value each and every customer, cherishing their diverse needs and taking immense pride in our capacity to serve a wide-ranging clientele. With over decades of unparalleled experience in the field of physiotherapy, we have honed a profound understanding of individual clinician nuances and the intricate demands of corporate entities.

Our suite of products and services is meticulously tailored to meet the precise requirements of each client, offering cutting-edge solutions for Ortho Rehab, Neuro Rehab, Cardio Rehab, and Sports Rehab. Whether you're initiating a new rehabilitation initiative or seeking to upgrade your current setup, our adept team is poised to provide expert guidance and support. Our unwavering commitment to exceeding customer expectations drives us to foster endur-

ing relationships and continually refine our offerings for optimal customer satisfaction.

Join us today and let us be your trusted partner in achieving excellence in the realm

of Physio-Rehab.

**Advanced Medical Equipment** 

# **1 MHz Heavy Duty Ultrasound** (Item Code: MDS – 004)

### **FEATURES:**

- Highly User-Friendly Control Panel
- Compact Design, Versatile Features
- Easy To Carry
- Good Quality Cover

### **ACCESSORIES:**

Ultra Sound Applicator :1

Ultra Sound Gel

: **1 Bottle** (250 ml)

**Carry Bag** 

:1

MEDICISATION DIGITAL ULTRASONIO MININA MININ	

<b>TECHNICAL SI</b>	PECIFICATION:
Input Supply	: 230 V Ac, 50 Hz
Size	: 8.75″ x 8″ x 3.5″ (LxWxH)
Weight	: 1.80 Kg. (Approx.)
Frequency	:1 MHz
Output mode	e : Continuous & Pulsed (1:1, 1:4, 1:7, 1:10)
Intensity	: Continuous Mode 0 to 3.5 watts/cm2 Pulsed Mode 1 .1 to 2.0 wattscm2
Timer	: 1 to 99 Minutes Programable
Buzzer	: Audio Buzzer at The End of Treatment Time

# **1 MHz Ultra Sound** (Ultra Z – Series) (Item Code: Z-501)

### **FEATURES:**

- Highly User-Friendly Control Panel
- •Wave Shape Biphasic
- Compact Design, Versatile Features
- Easy To Carry

### **ACCESSORIES:**

Carry Bag : 1

Ultrasound Gel : 1 Bottle (250 ML)

US Head & Cable :1



Input Supply	: 230 V Ac, 50 Hz
Size:	: 10.5" x 10.75" x 4.25" (LxWxH)
Weight :	: 1.80 Kg. (Approx)
Frequency	:1 MHz
Output Mode	: Continuous & Pulsed (1:1, 1:4, 1:7, 1:10)
Intensity	: Continuous Mode 0 to 3.5 watts/cm2
Pulsed Mode	: 1.1 to 2.0 watts /cm2
Timer	: 1 to 99 Minutes Programable
Buzzer	: Audio Buzzer at The End of Treatment Time

# **1 & 3 Mhz Ultrasound** (Ultra Z- Series) (Item Code: Z-502)

### **FEATURES:**

- 27 Pre-Programs Condition Wise
- Light LCD Display
- Manual Setting Option
- Highly User-Friendly Control Panel
- Compact Design, Versatile Features

### **ACCESSORIES:**

Ultra Sound Gel : 1 Bottle (250ML)

:1

Carry Bag



Input Supply	: 230 V Ac, 50 Hz
Size	: 10.5" x 12.5" x 4.75" (LxWxH)
Weight	: 1.80 kg (Approx)
Output	: Dual Head
Frequency	: 1 MHz & 3 MHz
Output Mode	: Continuous & Pulsed
Intensity	: 0.1 to 3.5 watt/cm2
Timer	: 1 to 99 Minutes Programable
Buzzer	: Audio Buzzer at The End of Treatment Time

# Four Channel TENS Auto (Item Code: MDS - 005)

### **FEATURES:**

- Highly User-Friendly Control Panel
- Universal & Economical Model
- Easy-to-Operate Design
- Six Types of Mode

### **ACCESSORIES:**

Electrode Cable	: 4
Rubber Electrode	:8
Velcro Set	:1
Ultra Sound Gel	: 1 Bottle (250 ml)
Carry Bag	:1



Input Supply	: 230 V Ac, 50 Hz
SIZE	: 8.5" x 10.75" x 3" (LxWxH)
Weight:	: 2.4 kg (Approx)
Output Voltage	e : 0-50 V with Load of 500 ohm
Pulse Width	: 50 to 350 micro sec
Frequency	: 02 to 150 Hz
Mode	: Six Types of Mode (Continuous, Brust1, Brust2, PWM, Sweep, Combi)
Timer	: 1 to 99 Minutes Programable
Buzzer	: Audio Buzzer at The End of Treatment Time

# **Four Channel TENS Auto** (Ultra Z- Series) (Item Code: Z-503)

### **FEATURES:**

- Highly User-Friendly Control Panel
- Universal & Economic Model
- Easy- To-Operate Design
- Six Types Of Mode

### **ACCESSORIES:**

Tens Wire	: 4
Round Rubber Electrode	: 8
Velcro Set	:1
Ultra Sound Gel	: 1 Bottle (250N
Carry Bag	:1



Input Supply	: 230 V Ac, 50 Hz
Size	: 10.5" x 10.75" x 4.25" (LxWxH)
Weight	: 1.85 kg. (Approx.)
Output Voltage	e : 0 to 50 V with Load of 500 ohm
Pulse Width	: 50 to 350 micro sec
Frequency	: 02 to 150 Hz
Mode	: Six Types of Mode (Continuous, Brust1, Brust2, PWM, Sweep, Combi)
Timer	: 1 to 99 Minutes Programable
Buzzer	: Audio Buzzer at The End of Treatment Time

# **5 IN 1 Combo** (125 Programs Computerized Machine) (Item Code: MDS-010)

### **FEATURES:**

- 125 Programs Condition Wise
- Backlight Lcd Display For IFT, TENS, & MS
- Separate Display For Timer, Intensity & Beat Frequency For Ultra Sound Therapy & Deep Heat Therapy
  - реер неат тнегару
- Simple & Beautiful Design

### **ACCESSORIES:**

Ultra Sound Applicator	:1
IFT+TENS+MS Wire	:4
IFT+MS+TENS Pad	:4
Deep Heat Pads	:2
Velcro Set : 1	:1
Ultra Sound Gel	:1 Bottle (250M)
Carry Bag	:1



### **TECHNICAL SPECIFICATION:**

Size	: 10" x 11.25" x 3.5" (LxWxH)
Weight	: 3.50 kg (Approx)
IFT Therapy mode	: 4P, 4PV, 2P
Carrier Frequency	: 2, 4KHz
Base	: 0-150 Hz
Sweep	: 0- 100Hz
Output current	: 0-100 mA
Ultrasound Therapy Modes	: Continuous and Pulse Ratios (1:1, 1:4, 1:7, 1:10)
Timer	: 1 TO 99 MINUTES Programmable
Buzzer	: Audio Buzzer at The End of Treatment Time

### **FEATURES:**

Interferential | Russian

Muscles Stimulator | TENS

Ultra Sound Therapy | Deep Heat

# Muscle Stimulator + Four Channel TENS Auto (Item Code: MDS-013)

### **FEATURES:**

- Highly User-Friendly Control Panel
- Compact Design, Versatile Features
- Easy To Carry
- Universal & Economical Unit

### **ACCESSORIES:**

TENS WIRE: 4

ROUND RUBBER ELECTRODE : 8

MS WIRES & PADS: 2 PC (WIRE), 2PC (PAD)

:1

:1

:1 Bottle (250ML)

Velcro Set /td> :1

Velcro Set

Ultra Sound Gel

Carry Bag



Input Supply	: 230 V Ac, 50 Hz
Size	: 8.5" x 10.75" x 3" (LxWxH)
Weight	: 2.40 kg (Approx., )
Output Voltage	e : 0 to 50 V with Load of 500 ohm
Pulse Width	: 50 to 350 micro sec
Frequency	: 02 to 150 Hz
Mode	: Six types of Mode (Continuous, Brust1, Brust2, PWM, Sweep, Combi)
Timer	: 1 to 99 minutes programmable
Buzzer	: Audio Buzzer at the end of treatment time

# **1 MHz Ultra Sound + Four Channel TENS** (Item Code: MDS-012)

### **FEATURES:**

- Independent US & FOUR Channel TENS
- Portable & Compact Model
- Soft Touch Buttons
- Adjustable Frequency
- SIX TYPES OF MODES FOR TENS
- 1 MHz ULTRASOUND
- Best For Different Chronic & Other Pain Management

### **ACCESSORIES:**

Tens Wires	: 4 Set
Round Rubber Electrodes	: 8 Pcs
US Head & Cable	:1
Velcro Strap	:1
Ultra Sound Gel	: 1 Bottle (250ml)
Carry Bag	:1



Input Supply	230 V Ac, 50 Hz	
Size	8.5″x 10.75″ x 3″ (LxWxH)	
Weight	2.5 kg. (Approx)	
For Ultra Sound		
Frequency	1 MHz	
Output Mode	Continuous & Pulsed	
Intensity	3.5 watt/cm2	
For FOUR CHANNEL TENS		
Output Voltage	0 to 50 V with Load of 500 ohm	
Pulse Width	50 to 350 micro sec	
Frequency	02 TO 150 Hz	
Mode	Six Types of Mode (Continuous, Brust1, Brust2, PWM, Sweep, Combi)	
Timer	1 to 99 Minutes Programmable	
Buzzer	Audio Buzzer at The End of Treatment Time	

# Interferential therapy Unit 70 Programs + Ultra Sound (Item Code: MDS-008)

### FEATURES:

- 70 Programs Condition Wise
- Russian Current
- Universal & Economical Unit
- Easy-To-Operate Design
- Highly User-Friendly Control Panel

### **ACCESSORIES:**

Electrode Cable	: 4
Size Rubber Pads	: 4
US Head & Cable	:1
Velcro Straps	:1
Ultrasound Gel	: 1 Bottle (250ml)
Carry Baa	:1



Input supply	230 V. Ac, 50 Hz	
Size	8.75"X 10.5" X 3.5" (LxWxH)	
Weight	2.8 kg. (approx.)	
For IFT		
Russian Current	2.5 K Hz Carrier Frequency	
Modes	4P, 4PV, 2P,	
Carrier Frequency	2K Hz, 4k Hz. Base: 0-150Hz, Sweep: 0-100Hz	
Output Current	0-100mA	
Timer	1 To 99 Minutes Programmable	
For Ultrasound		
Frequency	1 MHZ	
Output Mode	Continuous & Pulsed (1:1, 1:4, 1:7, 1:10)	
Intensity	Continuous Mode 0 to 3.5 watts/cm	
Pulsed Mode	1.1 TO 2.0 Watts / Cm	
Timer	1 to 99 Minutes Programmable	
Buzzer	Audio Buzzer at The End of Treatment Time	

# Interferential Therapy Unit with 70 PROGRAMS (Item Code: MDS - 007)

### **FEATURES:**

- 70 Programs Condition Wise
- Russian Current
- Universal & Economical Unit
- Easy-To- Operate Design
- Highly User-Friendly Control Panel

### **ACCESSORIES:**

Electrode Cable	: 4
Two Size Rubber Pads	: 4
Velcro Straps	:1
Ultra Sound Gel	:1 Bottle (250ml)
Carry Bag	:1

MICRO CONTROLLI BASED COMBISTIM	I.F.T 70 PROGRAM	
	888	

Input supply	: 230 v. Ac, 50 hz
Size	: 9.75" x 12.5" x 3.5" ( L x W x H)
Weight	: 2.5 kg. (approx.)
Russian Current	: 2.5 K Hz Carrier Frequency
Modes	: 4P, 4PV, 2P,
<b>Carrier Frequency</b>	: 2 K Hz, 4 K Hz. Base: 0-150 Hz, Sweep: 0-100 Hz,
Output Current	: 0-100 mA,
Timer	: 1 to 99 Minutes Programable
Buzzer	: Audio Buzzer at the End of Treatment Time

# Impact IFT (Item Code: MDS-006)

### **FEATURES:**

- Easy-To-Operate Design
- Easy To Carry
- Universal Unit
- User-Friendly

### **ACCESSORIES:**

Electrode Cable : 4	: 4
Two Size Rubber Pads	: 4
Velcro Straps	:1
Ultra Sound Gel	: 1 Bottle (250m
Carry bag	:1



Input Supply	: 230 V. Ac, 50 Hz
Size	: 8.5″x 10.75″x 3″ (LxWxH)
Weight	: 2 kg. (Approx.)
Modes	: 2P, 4P, 4PV, 4PV 45, 4PV 90
Medium Frequency	: 4 K Hz, 2 K Hz
Carrier Frequency	: 0-50 Hz, 50-80 Hz, 80-100 Hz, 100-120 Hz, 120-150 Hz
Output Current	: 0–100 mA,
Timer	: 1 TO 99 Minutes Programable
Buzzer	: Audio Buzzer at the End of Treatment Time

# Mini IFT + Ultra Sound (Item Code: MDS-006)

### **FEATURES:**

- Highly User-Friendly Control Panel
- Compact Design, Versatile Features
- Easy To Carry
- Good Quality Cover

### **ACCESSORIES:**

- Us Head & Cable :1
- Two Size Rubber Pads : 4
- Velcro Straps : 1

Ultra Sound Gel : 1 bottle (250ml)

:1

Carry Bag



Input supply	230 V. Ac, 50 Hz
Size	8.5"X 10.75" X 3.5" (LxWxH)
Weight	2.5 kg. (approx.)
For IFT	
Modes	2P, 4P, 4PV, 4PV 45, 4PV 90
Carrier Frequency	2K Hz, 4k Hz. Base: 0-150Hz, Sweep: 0-100Hz
Medium Frequency	4 K Hz, 2 K Hz
Carrier Frequency	0-50 Hz, 50-80 Hz, 80-100 Hz, 100-120 Hz, 120-150 Hz
Output Current : 0-100 mA,	
Timer	1 To 99 Minutes Programable
Buzzer	Audio Buzzer at the End of Treatment Time
For Ultrasound	
Frequency	1 MHZ
Output Mode	Continuous & Pulsed (1:1, 1:4, 1:7, 1:10)
Intensity	Continuous Mode 0 to 3.5 watts/cm
Pulsed Mode	1.1 TO 2.0 Watts / Cm
Timer	l to 99 Minutes Programmable
Buzzer	Audio Buzzer at The End of Treatment Time

# Laser Therapy Z- 500 (Ultra Z- Series) (Item Code: Z-500)

### **FEATURES:**

- Protocol 45 Programs
- Most Advance, Powerful, Versatile, User-Friendly
- Safety Lock To Prevent Accidental Exposure | Red Visible Pointed Probe
- Large Graphic Lcd Display
- Probe Red Visible Laser Output Wave Length: (610±30) Nm
- Pre-Program Saved With The Name Of Disease
- Auto Probe Detection
- Ir685nm-100 Mw Pointed Probe
- Display Energy Level Give Energy

### **ACCESSORIES:**

Laser Protective Goggles	:2
Good Quality Cover	:1
Spare Key for Safety Lock	:2

Ir685nm-100 mw Pointed Probe : 1

# C E MedsEasy

Input Supply	: 230 V Ac, 50 Hz
Size	: L 11″ * W 12.5″ * H 4.25″
Weight	: 2.5 kg (Approx)
Class: 3 B Laser probe	: 685nm - 100 mW pointed (module)
Wave Length	: 610±30nm
Treatment time	: 1 - 60 minutes
Pulse Frequency	: 9 steps 10 Hz to 5KHz Sweep : 10 Hz to 5 KHz
Pulse width	: 100 us nominal Technical Data
Mains voltage	: 110V/230V, 50 or 60 Hz

# Mini TENS (Item Code: MDS - 001)

### **FEATURES:**

- Easy-To-Operate Design
- Easy To Carry
- Universal Unit
- User-Friendly

### **ACCESSORIES:**

Electrode	Cables	:2
-----------	--------	----

- Rubber Pads : 4
- Velcro Straps :1
- Ultra Sound Gel : 1 Bottle (250ml)
- Carry Bag : 1



Input supply	: 230 v. Ac, 50 Hz
Size	: 5.5" x 5.5" x 2.75" (LxWxH)
Weight	: 1 kg. (Approx.)
Model	: Analog
Output	: 2 channel (4 pads)
Frequency	: 2-250 Hz
Intensity	: Adjustable

# Mini TENS + MS Combo (Item Code: MDS - 002)

### **FEATURES:**

- Dual Independent Channel TENS And Single Channel MS
- MS Currents: Galvanic & Faradic
- Portable & Compact Model
- Easy To Carry
- Light Weight

### **ACCESSORIES:**

Electrode Cables	:4
Square Rubber Electrodes	:4
Carry Bag	:1
Pen Electrode	:1
Velcro Set	:1
Ultra Sound Gel	:1Bc



Input Supply	: 230 V Ac, 50 Hz
Size	: 5.5" x 5.5" x 2.75" (LxWxH)
Weight	: 1 kg. (Approx.)
Model	: Analog
TENS Channel	: 2 Channel (4 Pads)
MS Channel	: Single Channel (2 Pads)
Treatment Mode (MS)	: Faradic & Galvanic
Timer	: 1 to 99 Minutes Programmable
Buzzer	: Audio Buzzer at The End of Treatment Time

# **Mini Ultrasound** (Item Code: MDS-003)

### **FEATURES:**

- Highly User-Friendly Control Panel
- Compact Design, Versatile Features
- Easy To Carry
- Good Quality Cover

### **ACCESSORIES:**

Ultra Sound Applicator : 1 Ultra Sound Gel : Bottle (250ml.) Carry Bag : 1



Input Supply	: 230 V Ac, 50 Hz
Size	: 5.5" x 5.5" x 2.75" (LxWxH)
Weight	: 1 kg. (Approx)
Frequency	: 1MHz
Output Mode	: Continuous & Pulsed
Intensity	: 3.4 watt / cm2
Timer	: 1 to 99 Minutes Programmable
Buzzer	: Audio Buzzer at The End of Treatment Time

# **Dual Channel Muscle Stimulator** (Item Code: MDS - 014)

### **FEATURES:**

- Highly User-Friendly Control Panel
- Compact Design, Versatile Features
- Easy To Carry
- Universal & Economical Unit

### **ACCESSORIES:**

Electrode Cables	: 4
Two Size Rubber Pads	: 4
Velcro Straps	:1
Ultra Sound Gel	: 1 Bottle (250ml
Carry Bag	:1



Input Supply	: 230 V Ac, 50 Hz
Size	: 8.5" x 10.75" x 3" (LxWxH)
Weight	: 2.40 kg. (Approx.)
Treatment Mode	: Faradic & Galvanic IG
Mode	: Child & Adult
Pulse Width	: 100 Milli Sec.
Frequency	: 1 Per Sec.
SF	: 1 TO 9 Sec. Duration
Timer	1 to 99 Minutes Programable
Buzzer	Audio Buzzer at The End of Treatment Time

# **Four Channel Digital Muscle Stimulator** (Item Code: MDS - 015)

### **FEATURES:**

- Highly User-Friendly Control Panel
- Compact Design, Versatile Features
- Easy To Carry
- Universal & Economical Unit

### **ACCESSORIES:**

Electrode	Cables	:4
-----------	--------	----

- Two Size Rubber Pads : 4
- **Velcro Straps** :1

**Carry bag** 

Ultrasound Gel : 1 Bottle (250ml)

:1



: 230 V Ac, 50 Hz
: 8.5" x 10.75" x 3" (LxWxH)
: 2.5 kg. (Approx.)
: Faradic & Galvanic IG
: 100 milli sec.
: 1 Per sec.
: 1 to 9 sec. Duration
: 1 to 99 Minutes Programable
: Audio Buzzer at The End of Treatment Time

# **Deep Heat** (Item Code : MDS - 011)

### **FEATURES:**

- It Is A Helpful In Relieving Pain
- Well As Treats Cellulite.
- 2 Channel Deep Heat Machine
- Used In Reduction Of Fat (Mainly Unwanted Cellulite) Deposits On Body
- Temperature Controlled Via Thermostat
- We Are Offering A Comprehensive Range Of Cellulite Deep Heat Therapy.
- This Deep Heat Therapy Is Designed Machine To Increase The Extensibility Of Soft Tissues, Remove Toxins From Cells, Enhance Blood Flow, Increase Functioning Of Tissue Cells, Encourage Muscle Relaxation And Help Relieve Pain
- Same Is Also Used In Reduction Of Fat (Mainly Unwanted Cellulite) Deposits On Body

### **ACCESSORIES:** 1

### **Carry Bag**



# Paraffin Wax Bath (Item Code: MDS – 016)

### **FEATURES:**

- Easy Mobility With Caster
- Highly Insulated Double Wall
- Automatic Temperature (30° To 110°) Control
- Available In Clinic And Hospital Model
- Outside Powder-Coated MS Tank



# Hydro Collator Pack (Item code : MDS - 017)

### **FEATURES:**

- Made Out Of Heavy Stainless-Steel Sheet
- Thermostat Control For Temperature Setting
- Highly Insulated Double Walled
- Four Moist Heat Pack Along With This.

### **ACCESSORIES:**

Four Moist Heat Pack



Input Supply	: 230 V Ac, 50 Hz
Power	: 2000 Watts
Watts Temperature	
Range	: 30 to 110
Weight	10KG (Approx.)

# Infra Ray Lamp (Item code : MDS - 024)

### **FEATURES:**

- Robustness
- Excellent Brightness
- Corrosion Resistance
- Easy To Use





# **New Product**



Ultrasound Gel (item code: 018)

SPECIFICATION: 250 ml | Water Soluble | Hypoallergenic | Salt Free





Hot & Cold Pack (Regular) (Item Code: MDS - 019)

Hot & Cold Pack (Large) (Item Code: MDS - 020)

**SPECIFICATION:** Long Lasting | Safe To Use | Eco Friendly

# **EXERCISE EQUIPMENTS**





**Assist Rolling Walkers** 

Join us today and let us be your trusted partner in achieving excellence in the realm of Physio-Rehab.



### **Contact Us**

- +91 97254 91989 | +91 90167 88893
- 🔀 Info@medseasy.in

www.medseasy.in

- 🕋 079-48975560
- FF-109-110, Seventh Heaven, 100, Ring Road, opp. Al Burooj, Makarba, Ahmedabad, Gujarat 380055

Note :- As our products are continuously enhanced through R & D, we reserve the right to change their specifications.

### An ISO 9001:2015, ISO 13485:2016, CE Certified Company